

Features

- Multi Assay Analyzer system
- Easy handling steps
- Short test time saves time and increases productivity
- Reliable, lab quality results



Afinion™ AS100 Analyzer

Simple ■ Fast ■ Reliable

- 1 Award winning design
 - 'Best Technology' at the 2011 European Mediscience Awards
 - Medical Design Excellence Awards (MDEA) competition. 2006
- 2 Intuitive operation and minimum training required
- 3 Fully enclosed testing environment increases user safety
- 4  USB connection for optional barcode reader and printer



Axis-Shield PoC Division
P.O. Box 6863 Rodeløkka
NO-0504 Oslo, Norway
www.axis-shield.com



Afinion™ AS100 Analyzer

Simple ■ Fast ■ Reliable

Connectivity

The Afinion™ Data Connectivity Converter (ADCC) transfers patient and control results from the Afinion™ AS100 Analyzer to LIS/HIS systems.

The results are transferred on the Ethernet in the standard formats HL7 or ASTM (selectable by configuration). The ADCC configuration is performed in a web interface that can be opened in a common web browser.

3-step procedure



Facts

- Size: 32 x 17 x 17 cm, weight: 5 kg
- LCD colour touch display
- No maintenance required
- No calibration needed
- Unique integrated failsafe system
- The Analyzer software can easily be upgraded via a USB flash drive to accommodate future tests

Available tests

Afinion™ HbA1c, Afinion™ ACR (albumin/creatinine ratio) and Afinion™ CRP

Please go to www.afinion.net to view a demo or to find your nearest distributor



Axis-Shield PoC Division
P.O. Box 6863 Rodeløkka
NO-0504 Oslo, Norway
www.axis-shield.com



afinion™