



# Sed-Chek<sup>®</sup> 2

The Industry Leading  
Erythrocyte Sedimentation Rate Control



Polymedco's Sed-Chek 2 is a whole blood reference control material designed to monitor ESR procedures. SED-CHEK 2 helps to monitor technique as well as environmental, physical and mechanical factors such as room temperature, tube position, vibration or instrument problems.

Sed-Chek 2 is an advanced formula and is the most comprehensive ESR control assayed for both automated and manual methods.

- ◆ Superior product performance
- ◆ Greatest stability of any ESR control - 18 months from date of manufacture
- ◆ 31 day open vial stability
- ◆ Room temperature storage
- ◆ Easy to mix formula
- ◆ Run like a patient sample - no voiding of diluent required
- ◆ Bi-level control that monitors drifts and trends in patient results due to environmental and procedural factors
- ◆ QC program with Peer Group Review available at no extra charge

## Sed-Chek 2 ESR Control

The number 1 selling ESR control in the world

The erythrocyte sedimentation rate test is a variable test where the result can be affected by environmental and procedural factors. Sed-Chek 2, a bi-level whole blood control, was designed for use with manual methods and automated ESR analyzers to monitor ESR procedures. Sed-Chek 2 employs a unique cell stabilizing process which mimics the rouleaux formation of patient samples.

Use of the Sed-Chek 2 ESR control will monitor technique as well as the environmental and physical factors which adversely affect ESR results. Improved laboratory confidence in accurate results is now assured by the use of the Sed-Chek 2 control.

Catalog Numbers	Description	Volume	Packaging
ESR-2CT	BI-LEVEL	48 mL	6 x 8 mL (3 x 8 mL each level)
ESR-2CTN	NORMAL	48 mL	6 x 8 mL
ESR-2CTA	ABNORMAL	48 mL	6 x 8 mL
For Sedimat 15			
ESR-5CT	BI-LEVEL	48 mL	6 x 8 mL (3 x 8 mL each level)
ESR-5CTN	NORMAL	48 mL	6 x 8 mL
ESR-5CTA	ABNORMAL	48 mL	6 x 8 mL
For StaRRsed <sup>®</sup> (Select Models)			
ESR-2SSN	NORMAL	30 mL	6 x 5 mL
ESR-2SSA	ABNORMAL	30 mL	6 x 5 mL

Please inquire about Polymedco's Sediplast<sup>®</sup> ESR System and Sedimat<sup>®</sup>15 automated ESR analyzer.

We offer flexible programs to make our total system available to your laboratory at affordable prices.

 **POLYMEDCO, Inc.**

510 Furnace Dock Road Cortlandt Manor, NY 10567

Phone: 800-431-2123 914-739-5400

Customer Service / Ordering Department Fax: 914-739-5890

Technical Service Fax: 914-739-0306

Reimbursement Support Hotline: Extension 228

[www.polymedco.com](http://www.polymedco.com)

[info@polymedco.com](mailto:info@polymedco.com)

**POLYMEDCO**

Sed-Chek 2 is a Registered Trademark of Polymedco, Inc.  
Sediplast and Sedimat are Registered Trademarks of LP ITALIANA  
StaRRsed is a Registered Trademark of R&R Mechatronics BV

© Copyright April 2002 Polymedco, Inc.  
PN 50179-00