

# QUICKVUE<sup>®</sup>

One-Step hCG  
Urine and Combo **TEST**

## Two Simple Tests You Can Depend On For The Early Detection of Pregnancy

### Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Positive may become visible sooner.

### Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 25 mIU/mL sensitivity in serum or urine.

### Convenient

- Room temperature storage.
- 18-month shelf life.

### Easy to Use


- Simple, one-step procedure uses only 3 drops of sample.
- Two-color endpoint for easy-to-read and interpret results.

### Reliable

- Highly sensitive assay with built-in controls.

### Lot-to-Lot Consistency

- Highly controlled and monitored manufacturing environment.



# QUIDEL<sup>®</sup>

For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at [www.quidel.com](http://www.quidel.com), or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

### Combo Procedure:

**Urine** 3

**Serum** 5

| Result | Urine     | Serum     |
|--------|-----------|-----------|
| +      | hCG COMBO | hCG COMBO |
| -      | hCG COMBO | hCG COMBO |
| ∅      | hCG COMBO | hCG COMBO |

### Urine Procedure:

3


| Result | Urine     |
|--------|-----------|
| +      | hCG URINE |
| -      | hCG URINE |
| ∅      | hCG URINE |



To order: Call your Quidel Account Manager at 1.800.874.1517.  
Outside the U.S. call 858.552.1100.

QuickVue hCG Urine – 25 Test Kit:  
Catalog # 20109  
hCG Control Set (Urine):  
Catalog # 00272

QuickVue hCG Combo – 50 Test Kit:  
Catalog #20110  
Serum hCG Control Set:  
Catalog #00281



For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at [www.quidel.com](http://www.quidel.com), or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

Distributor Information: