



Colo-Rectal- and **i**Colo-Rectal test. Easy, quick and sensitive.
Proven tests for the early recognition of bowel cancer.

- Two test methods for detecting occult blood in stools to ensure optimum results:
 - Colo-Rectal test:** Modified guaiac method
 - i*Colo-Rectal test:** Immunological test specific to human hemoglobin
(82% sensitivity and 99% specificity compared with colonoscopy)
- Hygienic, clean and quick
- Easy to use

Colo-Rectal- and **i**Colo-Rectal test. Step-by-step diagnostics for the early recognition of bowel cancer.

Colon cancer can be cured if it is detected early. Patients in stage T1 or T2 have the chance of surviving for 5 years. Unfortunately bowel cancer is often in its advanced stages before it begins to cause problems.

Tests for occult blood in stools are the only screening method whose efficiency has been proven in large-scale, randomized, controlled trials. By regularly repeating the test, the carcinoma can be detected at an earlier stage and the mortality rate can be reduced.

For this reason, Axonlab offers two test methods with high levels of sensitivity and specificity: the Colo-Rectal screening test and the immunological **i**Colo-Rectal test with diagnostic properties.



Colo-Rectal test

- Modified guaiac method
- Cost-effective, with high levels of sensitivity
- 3-day application
- Hygienic for patients and laboratory staff
- Recognized screening test
- Hygienic, clean and quick
- No special diet necessary, can be carried out during iron supplementation

Colo-Rectal test

0740284	5 tests
0740276	10 tests
0740268	50 tests


Stool collector option for Colo-Rectal tests

07033870	30 collectors
----------	---------------

NEW for both tests: Patient instructions in twelve languages – for increased ease of use.

German · French · Italian · English · Spanish
Portuguese · Dutch · Croatian · Slovenian
Polish · Hungarian · Turkish

Two test methods for optimum results




Usage and test evaluation: The Colo-Rectal and *i*Colo-Rectal tests are easy to use, even for patients. The tests can be evaluated after a short time by a doctor or a laboratory. The instructions are easy to understand and clearly laid out. Laboratory staff never come into direct contact with the stool sample.

*i*Colo-Rectal test

- Immunological test
- Specific detection of human hemoglobin
- High sensitivity of 82% and specificity of 99% compared with colonoscopy
- Hygienic and easy to carry out
- No direct contact with the stool sample

*i*Colo-Rectal test

07033310	10 tests
07033330	30 tests



Microbiology from Axonlab.

More than 300 media for medical and industrial microbiological diagnostics.

UrinAX TS and Urinax NOR

Bacterial count and resistance testing rolled into one.
Result obtained for resistance at the same time as the total bacterial count already after 18 to 24 hours.


AxonGram

Gram staining.
Stabilized color solutions,
long stability, individual
reagents available.




UrinAX

Immersion media.
Excellent growth and
easy bacterial count.



CountAX®

Hygiene control –
simply clean!
Total germ count (TPC)
and selective Entero-
bacteriaceae-Screen (E)
in surface fluids.



StorAX

Transport tubes
for stool samples.
Clean, odorless and
easy-to-use system for
treating parasites in stools.



ORIANA diag d.o.o.

www.oriana-diagnostics.com

Celjska cesta 15-17
SI-3252 Rogatec
Tel. +386 (0)3 818 61 12
Fax +386 (0)3 818 61 13
info@oriana-diag.com



Axon Lab AG

www.axonlab.com

Täferenstrasse 15
CH-5405 Baden
Tel. +41 (0)56 484 80 80
Fax +41 (0)56 484 80 99
info@axonlab.ch

Gewerbezone 1
A-6404 Polling
Tel. +43 (0)523 88 77 66
Fax +43 (0)523 88 77 66-12
info@axonlab.at